



HEALTHCARE
enhancing lifecare +



OPUS v12

The most powerful and versatile
monitoring with interactive interface



ENHANCING
LIFECARE +

KEY FEATURES

- + Standard: 5 leads ECG, Nellcor SpO2, Suntech NIBP, 2 TEMP, RESP, Li-ion battery, Thermal Printer
- + Optional: 3 leads ECG, Masimo SpO2, 4 IBP, EtCO2, C.O, AG, BIS
- + Interactive touchscreen up to 7 waveforms at a glance
- + 'Smart Sensor' automatically adjusting screen brightness
- + Rechargeable Li-ion battery with long life time
- + 120 full hours trend data storage
- + 48 hours of full disclosure waveforms review
- + OxyCRG, Drug Calculation
- + 26 types of Arrhythmia analysis, ST segment analysis and pacemaker detection
- + USB and SD memory data saving
- + Early warning alarm and customized alarm setting
- + Wire and wireless networking
- + External VGA monitor connection
- + Respiration and Phasein EtCO2 available
- + Rolling stands and wall mount available
- + Built-in integrated 3 channel thermal printer

Essential measurements deliver quick and accurate clinical data for enhanced patient care



'Smart Sensor' adjusts screen brightness automatically to maintain high legibility

Minimum seams of buttons make easy and effective disinfecting

User friendly interface with **12.1" touch screen**



KPI Healthcare Inc.

23865 Via del Rio, Yorba Linda California 92887 U.S.A.

Tel: +1 714 696 3434 | Email: info@kpimed.com | Web: www.kpihealthcare.com